


EN

REF Q-NCOV-01G

STANDARD Q COVID-19 Ag Test

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST


KIT CONTENTS



PREPARATION AND TEST PROCEDURE

■ PREPARATION

1. Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
2. Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
3. Check the test device and the desiccant pack in the foil pouch.



■ COLLECTION OF SPECIMEN

[Nasopharyngeal swab]

1. Insert a sterile swab into the nostril of the patient, swab over the surface of the posterior nasopharynx. Withdraw the sterile swab from the nasal cavity.



2. Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.




3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

4. Press the nozzle cap tightly onto the tube.

[Specimens in transport media]

1. Using a micropipette, collect the 350µl of specimen from the collection cup or VTM. Mix the specimen with an extraction buffer.
2. Press the nozzle cap tightly onto the tube.



[ANALYSIS OF SPECIMEN]

1. Apply 3 drops of extracted specimen to the specimen well of the test device.
2. Read the test result in 15-30 minutes.

Read in 15-30 mins.
Desiccant after 30 mins.

Desiccant after 30 mins.

■ LIMITATION OF TEST

1. The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.

2. The test should be used for the detection of SARS-CoV-2 antigen in human nasopharyngeal specimens only.

3. Neither the quantitative value nor the rate of SARS-CoV-2 antigen concentration can be determined by this qualitative test.

4. Failure to follow the test procedure and interpretation of results may adversely affect test performance and/or produce invalid results.

5. A negative test result may occur if the level of extracted antigen in a specimen is below the sensitivity of the test or if a poor quality specimen is obtained.

6. For more accuracy of immune status, additional follow-up testing using other laboratory methods is recommended.

7. The test result must always be evaluated with other data available to the physician.

8. A negative result may occur if the concentration of antigen or antibody in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or an molecular assay or ELISA.

9. Positive test results do not rule out co-infections with other pathogens.

10. Negative test results are not intended to rule in other coronavirus infection except the SARS-CoV.

11. Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.

■ Test principle

- STANDARD Q COVID-19 Ag Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and the test line are coated with anti-SARS-CoV-2 antibody.

- Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used as detectors for SARS-CoV-2 antigen. During the test, SARS-CoV-2 antigen in the specimen interact with monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-SARS-CoV-2 antibody. A colored band will be visible in the result window if SARS-CoV-2 antigen is present in the specimen. The intensity of colored test line will vary depending on the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents are working.

■ Kit contents

- ① Test device (individually in a foil pouch with desiccant)
- ② Extraction buffer tube
- ③ Nozzle cap
- ④ Sterile swab
- ⑤ Instructions for use

■ STORAGE AND STABILITY

- Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

■ WARNINGS AND PRECAUTIONS

1. Do not use the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use the extraction buffer tube of another lot.
4. Do not smoke, drink or eat while handling specimen.


* The presence of any line no matter how faint the result is considered positive.

* Positive results should be considered in conjunction with the clinical history and other data available.


Negative



Positive



Invalid



CAUTION

Do not read test results after 30 minutes. It may give false results.

INTERPRETATION OF TEST RESULT

* "C" Control Line | * "T" Test Line

Negative		
Positive		
Invalid		

1. A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).

2. A white band will appear in the bottom section of the result window. This band is test line of SARS-CoV-2 antigen (T).

3. Even if the control line is faint, or the test line is not uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.

* The presence of any line no matter how faint the result is considered positive.

* Positive results should be considered in conjunction with the clinical history and other data available.

■ KIT CASUTUSLAHUS

■ STANDARD Q COVID-19 Ag Test

■ STANDARD Q COVID-19 Ag Test